INSTRUCTION OF USE CORONAVIRUS TEST Novel Coronavirus (2019-nCoV) Antibody IgG/IgM Assay (Colloidal Gold)																								
INTENDED USE																								
The kit is used for the in vitro qualitative determination of Novel Coronavirus (2019-nCoV) IgG/IgM antibody in human serum, plasma or whole blood. It is used as an aid in the diagnosis of infection with Novel Coronavirus (2019-nCoV).It's not the only diagnostic criteria.																								
The important immunodiagnosis basis of Novel Coronavirus(2019-nCoV) infection are eliciting immune response and prompting B cell to create specific antibody. In earlier stage of infection, the detection aims at specific antibody IgM against nucleoprotein to find earlier infections. In infection and convalescent period, detecting antibody IgG. Continuously detection of antibodies and clinical observation which can be used to learn the regulation of virus elimination and large-scale epidemiological investigation. Therefore, the detection of Coronavirus test Novel coronavirus (2019-nCoV) antibodies IgG/IgM is a great significance to control the large-scale transmission of Novel Coronavirus (2019-nCoV) effectively																								
DESCRIPTION OF SYMBOLS USED																								
CORONAVIRUS TEST Novel Coronavirus (2019-nCoV) Antibody IgG/IgM Assay (Colloidal Gold) The following symbols may have been used in the labelling of this product.																								
<table> <tbody><tr> <td> EC</td> <td> REP</td> <td>Authorized Representative in Europe</td> <td> Manufacturer</td> </tr> <tr> <td> IVD</td> <td> LOT</td> <td>In Vitro Diagnostic Use</td> <td> Lot Number</td> </tr> <tr> <td> Expiry Date (Use by...)</td> <td> Temperature Limitation (store at 4~30°C)</td> <td></td> <td> Temperature Limitation (store at 4~30°C)</td> </tr> <tr> <td> Keep away from Sunlight</td> <td> Keep upright</td> <td></td> <td> Keep upright</td> </tr> <tr> <td> Do not reuse</td> <td> Read the instruction before proceeding</td> <td></td> <td> Read the instruction before proceeding</td> </tr> <tr> <td> Data of production</td> <td></td> <td></td> <td></td> </tr> </tbody></table>	 EC	 REP	Authorized Representative in Europe	 Manufacturer	 IVD	 LOT	In Vitro Diagnostic Use	 Lot Number	 Expiry Date (Use by...)	 Temperature Limitation (store at 4~30°C)		 Temperature Limitation (store at 4~30°C)	 Keep away from Sunlight	 Keep upright		 Keep upright	 Do not reuse	 Read the instruction before proceeding		 Read the instruction before proceeding	 Data of production			
 EC	 REP	Authorized Representative in Europe	 Manufacturer																					
 IVD	 LOT	In Vitro Diagnostic Use	 Lot Number																					
 Expiry Date (Use by...)	 Temperature Limitation (store at 4~30°C)		 Temperature Limitation (store at 4~30°C)																					
 Keep away from Sunlight	 Keep upright		 Keep upright																					
 Do not reuse	 Read the instruction before proceeding		 Read the instruction before proceeding																					
 Data of production																								

The kit is used for the in vitro qualitative determination of Novel Coronavirus (2019-nCoV) IgG/IgM antibody in human serum, plasma or whole blood. It is used as an aid in the diagnosis of infection with Novel Coronavirus (2019-nCoV).It's not the only diagnostic criteria.

The important immunodiagnosis basis of Novel Coronavirus(2019-nCoV) infection are eliciting immune response and prompting B cell to create specific antibody. In earlier stage of infection, the detection aims at specific antibody IgM against nucleoprotein to find earlier infections. In infection and convalescent period, detecting antibody IgG. Continuously detection of antibodies and clinical observation which can be used to learn the regulation of virus elimination and large-scale epidemiological investigation. Therefore, the detection of Coronavirus test Novel coronavirus (2019-nCoV) antibodies IgG/IgM is a great significance to control the large-scale transmission of Novel Coronavirus (2019-nCoV) effectively

DESCRIPTION OF SYMBOLS USED																								
CORONAVIRUS TEST Novel Coronavirus (2019-nCoV) Antibody IgG/IgM Assay (Colloidal Gold) The following symbols may have been used in the labelling of this product.																								
<table> <tbody><tr> <td> EC</td> <td> REP</td> <td>Authorized Representative in Europe</td> <td> Manufacturer</td> </tr> <tr> <td> IVD</td> <td> LOT</td> <td>In Vitro Diagnostic Use</td> <td> Lot Number</td> </tr> <tr> <td> Expiry Date (Use by...)</td> <td> Temperature Limitation (store at 4~30°C)</td> <td></td> <td> Temperature Limitation (store at 4~30°C)</td> </tr> <tr> <td> Keep away from Sunlight</td> <td> Keep upright</td> <td></td> <td> Keep upright</td> </tr> <tr> <td> Do not reuse</td> <td> Read the instruction before proceeding</td> <td></td> <td> Read the instruction before proceeding</td> </tr> <tr> <td> Data of production</td> <td></td> <td></td> <td></td> </tr> </tbody></table>	 EC	 REP	Authorized Representative in Europe	 Manufacturer	 IVD	 LOT	In Vitro Diagnostic Use	 Lot Number	 Expiry Date (Use by...)	 Temperature Limitation (store at 4~30°C)		 Temperature Limitation (store at 4~30°C)	 Keep away from Sunlight	 Keep upright		 Keep upright	 Do not reuse	 Read the instruction before proceeding		 Read the instruction before proceeding	 Data of production			
 EC	 REP	Authorized Representative in Europe	 Manufacturer																					
 IVD	 LOT	In Vitro Diagnostic Use	 Lot Number																					
 Expiry Date (Use by...)	 Temperature Limitation (store at 4~30°C)		 Temperature Limitation (store at 4~30°C)																					
 Keep away from Sunlight	 Keep upright		 Keep upright																					
 Do not reuse	 Read the instruction before proceeding		 Read the instruction before proceeding																					
 Data of production																								

CORONAVIRUS TEST Novel Coronavirus (2019-nCoV) Antibody IgG/IgM Assay (Colloidal Gold)
The following symbols may have been used in the labelling of this product.

 EC	 REP	Authorized Representative in Europe	 Manufacturer
 IVD	 LOT	In Vitro Diagnostic Use	 Lot Number
 Expiry Date (Use by...)	 Temperature Limitation (store at 4~30°C)		 Temperature Limitation (store at 4~30°C)
 Keep away from Sunlight	 Keep upright		 Keep upright
 Do not reuse	 Read the instruction before proceeding		 Read the instruction before proceeding
 Data of production			

PRINCIPLES OF THE PROCEDURE

CORONAVIRUS TEST Novel Coronavirus (2019-nCoV) Antibody IgG/IgM test is performed using the Novel Coronavirus (2019-nCoV) Antibody IgG/IgM Assay Kit on Colloidal gold-Immunochromatographic system. The detection card contained: 1) Recombinant antigen of Novel Coronavirus (2019-nCoV) labelled by Colloidal gold and Colloidal gold-labelled rabbit IgG antibody 2) Nitrocellulose membrane fixed with two detection lines(T1 line and T2 line) and one quality control line(C line).T1 line is fixed with monoclonal-anti human IgG antibody to test the antibody IgG of Novel Coronavirus (2019-nCoV). T2 line is fixed with monoclonal-anti human IgM antibody to test the antibody IgM of Novel Coronavirus (2019-nCoV).C line is fixed with goat anti rabbit IgG antibody.

Add defined amount of specimens to be tested into the sample hole of the detection card. The antibody IgG/IgM in specimens react with the Novel Coronavirus(2019-nCoV) antigen labelled by colloidal gold. After reacting, the immune complexes are formed. The complexes captured by the anti-human antibody IgG/IgM immobilized on the membrane, appearing a purple red T1/T2 line. This Show “Positive” for antibody IgG/IgM of Novel Coronavirus (2019-nCoV).

MAIN COMPONETS

- Detection card of Novel Coronavirus (2019-nCoV) antibody IgG/IgM 1 pc. Recombinant antigen of Novel Coronavirus (2019-nCoV) labeled by Colloidal gold, rabbit IgG beled by Colloidal gold, anti-human IgM, anti-human IgG, goat anti-rabbit IgG.
- Phosphate buffer solution. 1bottle × 1ml
- Dropper 1 pc
- Lancet 1 pc
- Alcohol pad 1 pc

STORAGE CONDITIONS AND EXPIRY DATE

The kit should be stored at 4°C~30°C and protected from light. The shelf-life of kit under these storage condition is 24 month. The detection card should be used within 1 hour after being opened. If the temperature is higher than 30°C or in high humidity (more than 60%) environment, it should be out of the box and wait before use until it reaches room temperature and humidity (15-30°C). See manufacturing date and expiry date on label.


SPECIMEN COLLECTION AND PREPARATION

- Serum, plasma, and whole blood samples can be used.
- There is no special requirement for the patient/donor before specimens collection, and fasting is not required.
- Collect specimens in accordance with standardized blood collection techniques, and the pre-treatments are conducted in accordance with standard procedures.
- Do not use any haemolyzed blood and chyle specimens to avoid misjudgement.
- Samples that have been contaminated with microorganisms cannot be used for testing.
- Avoid multiple freeze-thaw procedures. Specimens may be stored at 2°C -8°C for up to 3 days. For long-term storage, specimens should be kept below -20°C. Do not freeze the specimens of whole blood. If specimens cannot be detected without delay, they shall be kept at 2°C -8°C. Bring specimens to room temperature prior to testing. Frozen specimens must be mixed thoroughly after thawing before testing.

TESTING PROCEDURE

- Serum, plasma, and whole blood samples: Frozen specimens must be taken to room temperature (15°C-30°C) prior to testing, completely thawed and mixed well before testing.

The fingerstick collection of blood: Wipe the inside of the ring finger with an alcohol cotton swab to disinfect, position the blood sample device firmly against puncture site. . Hold safety lancet between fingers. Do not remove the device from the site until an audible click is heard. Position the blood sample device firmly against puncture site. To activate, press the blood sample device firmly against the puncture site.
- Remove the detection card and the dropper from the pouch and put it flat on the table.
- Label the detection card with specimens number.
- Hold pipette or dropper vertically and transfer specimens 10µl (1 drop with the dropper) to the sample hole on detection card. Then immediately add 2 drops with the bottle (approximately 50µl~70µl) of sample dilution. Ensure that no bubbles are generated during operation.
- Start the timer and read the results in 15 minutes.

Note: Do not take more than 15 minutes to read the results!


After observing and recording the results, please discard the detection card to avoid confusing the result judgment. Please take a picture if you want to save the result for a long time.

INTERPRETATION OF RESULTS


1. **Negative:** Only the control band is visible, and no band is visible on detection line T1, T2, showing no antibody of Novel Coronavirus (2019-nCoV) was detected. The result is negative (as shown below).


2. **Positive:**
2.1. The control band and T1 band are both visible, showing antibody IgG of Novel Coronavirus (2019-nCoV) was detected. The result is positive of IgG(as shown below).


2.2. The control band and T1 band are both visible, showing antibody IgM of Novel Coronavirus (2019-nCoV) was detected. The result is positive of IgM (as shown below).

2.3. The control band, T1 band and T2 band are all visible, showing antibody IgG/IgM of Novel Coronavirus (2019-nCoV) was detected. The result is positive of IgG/IgM(as shown below).

3. **Invalid**
If the Control band is invisible, whether the T1/T2 band is visible or not, the assay is invalid (as shown below). Repeat the assay with a new card.

Int	1. Although the test result shows "Negative", related symptoms remains. It is recommended to make other medical checks and diagnose the causes.
	2. The test result shows "Positive", further review is recommended immediately.
	3. The kit is suitable for preliminary screening of suspected cases. The final result must be interpreted together with doctor and other clinical information available to the physician.

LIMITATION OF TEST METHOD

- The test results obtained by other methods are not directly comparable with the measurement results of this product.
- Inaccurate results may occur due to operations and specimens. Repeat testing to ensure the consistent results.
- The kit is limited to the qualitative test and assisted diagnosis. It cannot detect the accurate concentration antibody of specimens.
- A negative result can occur during early stage of infection if the quantity of the antibodies IgG/IgM to Novel Coronavirus (2019-nCoV) in the specimen is below the detection limits or not develop antibodies. It should be go to hospital for examination. During the re-examination, the last specimen collected should be detected again to confirm whether there is a serological positive or a significant increase in titer.

PERFORMANCE CHARACTERISTICS

- Appearance: all components of the kit should be complete; labels shall be clear and easily legible; The liquid is no leakage and clear.
- Testing specimens with a titer of 1:320 in Novel Coronavirus (2019-nCoV) antibody IgG and IgM, no hook effect was observed.
- After preliminary evaluation, it was basically confirmed that the clinical performance of the kits could meet the outbreaks of emergency. More clinical will be collected after kit launch to confirm its clinical performance.
- The positive reference coincident rate of antibody IgG: positive reference coincident rate of enterprise is 10/10; The positive reference coincident rate of antibody IgM: positive reference coincident rate of enterprise is 10/10.
- The negative reference coincident rate: positive reference coincident rate of enterprise is 10/10.
- The limit of detection (LOD): For IgG antibody detection, S6 should be negative, no requirements with S5; S1, S2, S3,S4 should be positive; For IgM antibody detection, L6 should be negative, no requirements with L5, L1, L2, L3,L4 should be positive.
- Repeatability

- 6.1) Taking two repetitive references of antibody IgM for testing. Testing each reference for 10 times. All results should be positive and uniform in color.
- 6.2) Taking two repetitive references of antibody IgM for testing. Testing each reference for 10 times. All results should be positive and uniform in color.
7. Analytical Specificity
 - 7.1) Cross-reaction: The antibody of Novel Coronavirus HKU1, antibody of Novel Coronavirus OC43, antibody of Novel Coronavirus NL63, antibody of Novel Coronavirus 229E, antibody of influenza virus, antibody of respiratory syncytial virus(RSV), antibody of adenoviridae, antibody of enterovirus had no Cross-reaction with this kit.
 - 7.2) Interfering substances: The interfering substances such as bilirubin, haemoglobin, triglyceride, rheumatoid factor, etc. showing neither inhibition nor enhancements.

PRECAUTIONS

- The kit is one-time use in vitro detection kit, please read the INSTRUCTION of use carefully before testing. Do not reuse!
- If the aluminum foil bag of the test card is damaged, it should be discarded.
- Ensure that a proper amount of specimens are used for testing. Too much or too little specimens may lead to deviations in results.
- The detection cards stored at low temperature, bring them to room temperature (15°C ~ 30°C) prior to opening and testing. To avoid moisture absorption.

REFERENCES

- National Health and Health Commission “Diagnosis and Treatment of Pneumonia of Novel Coronavirus Infection” (Trial version 3)
- World Health Organization: Clinical management of severe acute respiratory infection when Novel Coronavirus (nCoV) infection is suspected: Interim Guidance.

 CE	 REP	SUNGO Europe B.V. OlympischStadion 24, 1076DE Amsterdam, Netherlands Tel: +31 20 21 11 106 Fax: +31 20 21 11 106	 Importer/Dystrybutor ESJOT Sp. z o.o. Graniczna 2, 54-610 Wrocław, Poland info@coronavirus-test.com.pl
 Manufacturer			

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